

RESEARCH ARTICLE

Polymer-spray mass spectrometric detection and quantitation of hydrophilic compounds and some narcotics

Maria T. Dulay  | Richard N. Zare 

Department of Chemistry, Stanford University, Stanford, CA 94305, USA

Correspondence

R. N. Zare, Department of Chemistry, Stanford University, Stanford, CA 94305, USA.
Email: zare@stanford.edu

Funding information

National Institute of General Medical Sciences, Grant/Award Number: 1R41 GM113337-01; National Institute of Mental Health, Grant/Award Number: R01-MH112188

Rationale: High-throughput screening of biofluids is essential in monitoring concentration of a variety of drugs to determine their efficacy and toxicity. Organosiloxane polymers prepared by sol-gel chemistry as sample supports, and electrospray ionization emitters in a single material and as an alternative to paper substrates, is described in this study.

Methods: Hydrophobic drugs and hydrophilic streptomycin were analyzed by polymer-spray mass spectrometry with an LTQ-Orbitrap mass spectrometer. Drug samples in urine (1–2 μL) were deposited on an OSX polymer, allowed to dry, then electrosprayed from the polymer tip into the mass spectrometer without sample pretreatment. The OSX polymers, whose polarity and porosity can be controlled, were prepared by sol-gel chemistry where methyl-substituted alkoxy silanes were hydrolyzed in the presence of a pore template and an acid catalyst.

Results: Five nanograms each of seven narcotic drugs were detected in <1 min (relative standard deviation (RSD) of response <1% for each drug). Calibration curves of cocaine and streptomycin in urine were used to establish the performance of the polymer. For sample 1 ($n = 2$), the mean recovery for cocaine was 81% with paper and 90% with polymer. Streptomycin is detected with polymer, not with paper; for samples 1 and 2 ($n = 3$), mean recovery was 97% and 95%, respectively.

Conclusions: Organosiloxane polymers achieve more sensitive analysis than paper, allowing for more accurate quantitation of both hydrophobic and hydrophilic drug compounds. The ability to tailor the polymer polarity and porosity allows for the synthesis of a wide range of polymers, and thus opens many possibilities for further development and applications.

1 | INTRODUCTION

Improving the management of a variety of diseases caused by pathogenic organisms relies on the correct dosage of antibiotic drugs for a patient. Measuring the concentration of antibiotics administered to a patient as part of the dosing regimen is often done by repeated measurements of the drug concentration over time to determine the therapeutic range to sustain drug efficacy, avoid drug-related toxicity, and decrease multiple-drug resistance. Streptomycin, the most widely used aminoglycoside antibiotic used to treat infections in humans, is a highly polar (hydrophilic) compound with narrow therapeutic index (20–30 $\mu\text{g}/\text{mL}$) and toxicity (50 $\mu\text{g}/\text{mL}$), which may be severe or irreversible.¹ Between 29% and 89% of the administered dose is excreted unchanged in the urine within 24 h.² Streptomycin is also used as a growth promoter in food-producing animals with the potential to develop streptomycin-resistant pathogens.³ The detection and quantitation of streptomycin in biofluids often requires the use of

hydrophilic interaction chromatography (HILIC), an alternative form of high-performance liquid chromatography (HPLC), where uncharged polar compounds are separated on polar stationary phases coupled with mass spectrometry⁴ or matrix-assisted laser desorption/ionization time-of-flight mass spectrometry (MALDI-TOF MS).⁵ However, these analytical methods suffer from increased time of analysis because of sample pretreatment and labor-intensive chromatographic procedures. Rapid analytical methods based on MS for high-throughput analyses would provide important information on proper dosing of these drugs during antibiotic therapy and their detection in commercial food products.

Electrospray ionization mass spectrometry (ESI-MS) has been shown to be an invaluable tool in the analysis of chemical compounds in a variety of different matrices, including biological fluids. In conventional ESI, a small-diameter capillary is used to generate the electrospray that results in the formation of analyte ions. However, one of the drawbacks of using a capillary is the danger of clogging.

Over the past two decades, improvements to sampling and ionization have been made, including the development of ambient ionization techniques and ESI on solid substrates, such as metals and paper.⁶ These improvements have led the way for the development of paper-spray mass spectrometry, which is an appealing form of ambient ionization mass spectrometry because only a small volume of analyte, typically some biofluid, is deposited on the surface of the paper support, allowed to dry, and then electrosprayed from the tip of the paper into a mass spectrometer.⁷⁻¹⁰ This technique, which requires no pneumatic assistance to transport analytes to the mass spectrometer inlet, can be beneficial for use in on-site therapeutic drug monitoring in biofluids and in monitoring the quality of food products and can allow for high-throughput analyses. Open-surface sampling avoids the clogging problem of capillaries and allows for more choices of sampling support materials, such as a synthetic polymer.

An important feature of sampling materials for ambient ionization mass spectrometry is their polarity. Hydrophobic materials such as polymethylmethacrylate and polytetrafluoroethylene have been shown to have enhanced signal intensity and stability over less hydrophobic materials used in desorption ionization techniques.^{11,12} Paper coated with a variety of different compounds has also been shown to have enhanced performance as sample support in paper-spray mass spectrometry.¹³⁻²² We have found that hydrophobic polymer sample supports allow for the detection of hydrophilic drugs, such as streptomycin, that are not easily detected, if at all, using paper-spray mass spectrometry.

We describe here a modification of paper-spray mass spectrometry in which the paper support is replaced by a suitable organic-inorganic hybrid organosiloxane (OSX) polymer whose porosity and chemical composition can be carefully controlled by use of a sol-gel process,²⁴ allowing for facile modification of the polymer surface (e.g., control of the surface hydrophobicity). The sol-gel process offers a simple, easy, and inexpensive route to make these OSX polymers. OSX polymers are used for the first time as a sampling support and an ESI emitter in a single material for the detection and quantitation of hydrophobic narcotics ($\log P \sim 0.9$ to 2.3, where $\log P$ is the partition coefficient of the concentration of the compound in octanol to water),²⁴ and for the detection and quantitation of streptomycin ($\log P -6.4$),²⁴ a hydrophilic antibiotic drug, which is not detectable using paper-spray mass spectrometry.

2 | EXPERIMENTAL

2.1 | Chemicals

Methyltrimethoxysilane (MTMS), streptomycin, dihydrostreptomycin, lidocaine, narcotic drugs and their deuterated analogs were purchased from Sigma-Aldrich (St Louis, MO, USA) and used without further purification. Narcotic drugs were purchased as 1 mg/mL in methanol or acetonitrile: amphetamine, methamphetamine, 3,4-methylenedioxyamphetamine (MDA), 3,4-methylenedioxymethamphetamine (MDMA), 3,4-methylenedioxyethylamphetamine (MDEA), morphine, and cocaine. All solvents (MS-grade) and hydrochloric acid were purchased from Thermo Fisher Scientific (Waltham, MA, USA).

2.2 | Preparation of organosiloxane (OSX) polymer

OSX polymer was prepared by stirring 2.5 mL methyltrimethoxysilane (MTMS) with 0.72 mL water, 0.18 mL toluene, and 0.035 mL concentrated HCl for 15 min at room temperature. The reaction solution was cast into the wells of a 24-well polystyrene multiculture plate at various volumes, ranging from 0.75 mL to 1.5 mL, then aged and cured at room temperature for 21 h. The polymers were cut into triangles (base ~ 1.3 cm, each of the two sides ~ 1.6 cm, tip angle $\sim 49^\circ$) using conventional scissors. The polymer product was cleaned of any unreacted starting materials by immersion in acetonitrile for 2 h followed by immersion in water for approximately 18 h. The polymer triangles were dried under ambient conditions before use.

2.3 | Sample preparation

Each drug stock solution was prepared by dissolving the drug in water. Drug-spiked urine samples were prepared by directly adding the appropriate volume of drug stock solution into urine.


2.4 | Calibration curves

Standard samples of streptomycin were prepared in urine. Linear regression analysis was carried out on known added concentrations of streptomycin against the peak intensity of streptomycin to dihydrostreptomycin (internal standard). Standard samples of cocaine in urine were prepared in a similar manner, and linear regression analysis was carried out using isotope-labeled cocaine-d3 as the internal standard. The regression coefficient (R^2), slope, and intercept, and equation of the resulting standard curves were determined. Calibration curves for streptomycin in urine for polymer spray, and cocaine in urine for both paper and polymer spray, are presented in the supporting information.

2.5 | Polymer-spray mass spectrometry

Scheme 1 is a diagram of the polymer-spray setup. Sample volumes of 1 to 2 μL were loaded near the tip of a polymer or paper triangle. The sample droplet was dried under ambient conditions for 1 h. The tip of the polymer triangle was placed ~ 10 mm from the mass spectrometer inlet. All polymer materials used in the experiments were prepared by casting 0.75, 1.0, or 1.25 mL of reaction volume into molds. High voltage (5 kV) was applied by connecting the back end of the triangle to the HV power supply of the mass spectrometer via a flat stainless steel alligator clip. For paper-spray experiments, Whatman filter paper (Grade 1; Fisher Scientific, Waltham, MA, USA) was used and the same process used for polymer-spray was used for paper-spray. The spray solvent was a 1:1 (v/v) or 9:1 (v/v) combination of methanol and water with 0.1% by volume of formic acid.

Mass spectrometry was performed with an LTQ-Orbitrap XL (Thermo Fisher Scientific, Waltham, MA, USA). Mass spectra were acquired across m/z 100–500 for narcotic drugs and m/z 200–640 for streptomycin. The capillary temperature was set to 275°C, and the spray voltage was 5 kV. All mass spectra were acquired in positive-ion mode.



SCHEME 1 Diagram of polymer-spray mass spectrometry

Tandem mass spectrometry (MS/MS) of streptomycin was performed on the molecular ion m/z 582.27 ($[M + H]^+$) for structural confirmation through collision-induced dissociation (CID) with helium as the collision gas in the CID cell. The normalized collision energies were set to between 15% and 30% in the CID cell.

3 | RESULTS AND DISCUSSION

3.1 | OSX polymer

A methyl-substituted alkoxy silane reagent was specifically chosen to prepare the OSX polymers because the presence of methyl groups in the polymer network produced a mechanically robust polymer that avoids shrinkage and cracking during the ambient evaporation of solvent from the polymer pores. Nakanishi and co-workers²⁵ reported on the role that methyl groups in a sol-gel polymer play in creating a "spring back" effect unlike with polymers made with tetraalkoxysilanes that are commonly observed to crack during ambient drying. Toluene served as a pore template and as a solvent to control phase separation of the polymer precursors during the hydrolysis and condensation steps of the sol-gel process. Other pore templates, such as polyethylene glycol and CTAC, used in the preparation of OSX polymers were strongly imbedded in the polymer network and could not be completely removed during the cleanup procedure. These templates strongly interfered with the mass spectrometric detection of the drug compounds and therefore they were avoided as pore templates. The hydrophobic OSX polymer used in these experiments had a water contact angle of $\sim 95^\circ$. A drop of pure methanol deposited on the surface spread across the polymer. Because of the hydrophobic nature of the OSX polymer, spreading of an aqueous solution is minimized as compared to paper, which is hydrophilic. Figure 1 shows a comparison of the spreading of a 1- μ L droplet of methylene blue in water deposited on the surface of filter paper and OSX polymer. The methylene blue-water droplet is absorbed by the filter paper, resulting in uniform distribution of the dye at the top and the back sides of the filter paper. However, a similar volume appears as a small droplet on the OSX top surface, maintaining a higher surface tension in response to the hydrophobic surface of the polymer, and thereby concentrating the sample to a smaller surface area than on filter paper. Methylene blue is not observed on the back side of the OSX polymer, indicating that the dye did not penetrate through the polymer.




FIGURE 1 Photos of 1- μ L droplet of methylene blue in water on an OSX polymer (left) and filter paper (right) taken just after deposition of the droplets [Color figure can be viewed at wileyonlinelibrary.com]

The thickness of the polymer increased with increasing volume of reaction solution used in casting of the polymer in the molds (Table 1). However, the thickness of the polymer did not appear to affect the performance of the polymer, although more care was needed to align the thicker polymers to the inlet of the mass spectrometer.

The polymers were characterized by scanning electron microscopy (SEM). In Figure 2A, the surface topography of the polymer is comprised of circular divots of diameters up to 10 μ m. Unlike porous filter paper having an average pore size diameter of 11 μ m, the exterior of the OSX polymer appears to have little macroscopic porosity. The polymer network, however, contains a range of pores with diameters up to 10 μ m, as seen in Figure 2B. The connectivity of the pores is

TABLE 1 Average thickness of OSX polymers as a function of volume of reaction solution casted

Sample support	Volume (mL)	Avg thickness (mm)
OSX	0.75	0.2 (n = 1)
OSX	1.0	0.42 \pm 0.06 (n = 9)
OSX	1.25	0.58 \pm 0.07 (n = 5)
Paper	-	0.1 \pm 0.02 (n = 3)




FIGURE 2 Scanning electron micrograph of OSX polymer, showing A, the surface topography and B, the pores within the polymer network

unknown as it could not be determined from the SEM images. Nitrogen adsorption porosimetry measurements of OSX polymers could not provide surface area or pore volume information, indicating that the pores are very tiny and not interconnected.

Signal was produced only when the polymer was wetted with the spray solvent and voltage was applied to the polymer. A 20- μ L volume of aqueous methanol spray solvent remained stationary on an OSX polymer surface in the absence of an applied spray voltage (Figure 3A). When a 5-kV voltage was applied to the wetted OSX polymer, a stable spray mist was observed from the tip of the polymer triangle, as recorded by an optical camera (Figure 3B), similar to what has been observed in droplet spray ionization using a glass microscope coverslip where the spray solvent elutes along the top surface

of the glass.²⁶ When compared to paper under the same spray conditions, the OSX polymer exhibited significantly increased signal stability. Figure 4 shows the signal of an OSX polymer lasting about 6.5 min as compared to paper where the signal decreased around 0.3 min.

3.2 | OSX polymer-spray of hydrophobic narcotic drugs

Figure 5 demonstrates the utility of polymer-spray MS for the detection of seven narcotic drugs (1–7) and their deuterated analogs (a–g) in urine, which is achieved in less than 1 min for 5 ng loading of each drug and deuterated drug. The repeatability of polymer-spray for each narcotic drug using a hydrophobic OSX polymer triangle for




FIGURE 3 A, OSX polymer with 20 μ L of 1:1 (v/v) methanol/water spray solvent deposited on its surface and aligned in front of the inlet of the mass spectrometer without any applied spray voltage. B, a spray mist from an OSX polymer with 20 μ L of a 9:1 (v/v) methanol/water spray solvent on the polymer surface with an applied voltage of 5 kV [Color figure can be viewed at wileyonlinelibrary.com]




FIGURE 4 TIC showing signal stability using 20 μ L spray solvent, 1:1 (v/v) methanol/water with 0.1% (v/v) formic acid for A, OSX polymer and B, paper. Spray voltage: 5 kV




FIGURE 5 Mass spectra of narcotic drugs (1–7) and their deuterated analogs (a–g) in urine detected by polymer spray. 9:1 (v/v) methanol/water spray solvent; 5 kV spray voltage; 1 μ L sample loading. Drug compounds: 1, amphetamine; 2, methamphetamine; 3, MDA; 4, MDMA; 5, MDEA; 6, morphine; 7, cocaine; a, amphetamine-d6; b, methamphetamine-d5; c, MDA-d5; d, MDMA-d5; e, MDEA-d6; f, morphine-d6; and g, cocaine-d3. Background peaks are marked by B

each run is shown in Table 2. The response of each drug is measured as the peak intensity ratio of a drug to its deuterated analog. The RSD of response for each drug is less than 1%.

Two different samples of cocaine in urine were analyzed and quantified using paper- and polymer-spray MS. The calibration curves obtained for both paper- and polymer-spray covered a linearity range of 0.033 to 33 μ g/mL of cocaine in urine. The standard curves demonstrated regression correlation coefficients (R^2) of 0.999 for both paper and polymer. Analyses of cocaine in samples 1 and 2 are shown in Table 3. The analysis sensitivity of cocaine with the OSX polymer was approximately 10% better than that for paper-spray as observed in quantitation experiments and the peak intensities higher by a factor

TABLE 2 Response (peak intensity ratio of drug:Deuterated drug) for seven narcotic drugs using an OSX polymer for each run ($n = 8$) electro sprayed at 5 kV with 9:1 (v/v) methanol/water spray solvent with 1- μ L sample loading of 5000 ng/mL concentration of each drug in urine

Drug	Avg Response	RSD, %
Amphetamine	0.907	0.027
Methamphetamine	1.05	0.02
MDA	1.08	0.01
MDMA	1.14	0.16
MDEA	1.23	0.01
Morphine	1.17	0.23
Cocaine	0.971	0.006

TABLE 3 Quantitation of two samples of cocaine in urine by paper- and polymer-spray MS electro sprayed at 5 kV with 9:1 (v/v) methanol/water spray solvent

Substrate	Sample 1 ($n = 2$) μ g/mL	Sample 2 ($n = 3$) μ g/mL
Paper	13.2 \pm 0.9	4.16 \pm 1.22
OSX polymer	14.7 \pm 1.6	4.07 \pm 0.11
Actual concentration	16.3	3.67

of 2 compared to paper-spray. Mean recoveries of 90.2% (RSD 10.9%) and 81.0% (RSD 6.8%) were obtained for sample 1 by polymer- and paper-spray, respectively. The calculated concentrations for cocaine in sample 2 by polymer- and paper-spray were higher in both cases, but the %RSD values were 2.7% and 29%, respectively. The favorable elution of the drugs from the OSX polymer surface, and therefore the higher analysis sensitivity of polymer-spray for the quantitation of cocaine in urine, may be attributed to (1) the localization of the cocaine sample to a smaller area of the polymer surface thereby concentrating the sample as compared to paper where absorption of the sample results in spreading across a wider area, as seen in Figure 1 for methylene blue in water, and (2) analytes may be 'trapped' within the entangled cellulose fibers that comprise the highly porous paper substrate, preventing all the analytes being sprayed into the mass spectrometer. When paper was coated with silica, improvements in the detection of pesticides in milk by paper-spray MS were reported.¹⁴

3.3 | Qualitative analysis of streptomycin by OSX polymer-spray

Hydrophobic drug compounds, such as the narcotic drugs described above, have positive logP values and are detected by both paper- and polymer-spray MS. In contrast, hydrophilic ones, such as streptomycin, have negative logP values. The logP value of streptomycin is -6.4 .²⁴ Unlike the narcotics, the hydrophilic streptomycin is not detected using paper-spray MS, but is detected by polymer-spray MS under the same conditions. With cellulose being a polysaccharide, its surface is dense with hydroxyl groups. It is believed that the polar groups of streptomycin interact strongly with the hydroxyl groups of the cellulose fibers, limiting their elution from the paper during the electrospray process. The hydrophobic methyl-containing OSX polymer may contain less hydroxyl groups, making for weaker affinity of streptomycin to the polymer surface, which may contribute to the more favorable elution efficiency of streptomycin in polymer-spray MS. Figure 6 shows the total-ion current (TIC) and mass spectra of 1000 μ g/mL streptomycin deposited at a volume of 2 μ L (2- μ g loading). As seen in Figure 6, streptomycin peaks are observed when a spray voltage of 5 kV is applied to the OSX polymer wetted with 20 μ L of 1:1 (v/v) methanol/water spray solvent: the singly and doubly charged ions, $[M + 2H]^{2+}$ (m/z 291.63) and $[M + H]^+$ (m/z 582.27), arise from protonation and the singly and doubly charged ions are from the hydrated form of the aldehyde group $[M + H_2O + H]^+$ (m/z 600.28) and $[M + H_2O + 2H]^{2+}$ (m/z 300.64), respectively.²⁷ The methanol adduct, $[M + MeOH + H]^+$ (m/z 614.29), is also observed.

Figure 7 shows full mass spectra obtained by depositing 1 or 2 μ L of a sample mixture of streptomycin (50 μ g/mL) and the hydrophobic drug lidocaine (50 μ g/mL, logP 2.44)²³ for comparison onto paper (Figure 7A) and OSX polymer (Figures 7B and 7C) and applying 5 kV spray voltage and 20 μ L 1:1 (v/v) methanol/water spray solvent. In Figure 7A, the $[M + H]^+$ peak at m/z 235.2 for lidocaine is observed, but no peaks can be attributed to streptomycin at a low sample loading of 50 ng. When the drug sample is electro sprayed using the OSX polymer, however, the streptomycin peaks are easily




FIGURE 6 TIC and mass spectrum of streptomycin (1000 $\mu\text{g}/\text{mL}$) in water deposited on OSX polymer. Spray solvent: 1:1 (v/v) methanol/water with 0.1% (v/v) formic acid. Spray voltage: 5 kV. Sample loading: 2 μL




FIGURE 7 Mass spectra of a mixture of lidocaine (50 $\mu\text{g}/\text{mL}$) and streptomycin (50 $\mu\text{g}/\text{mL}$) in water deposited on A, paper, 1 μL loading; B, OSX polymer, 2 μL loading; and C, OSX polymer, 1 μL loading. Lidocaine: 1 $[\text{M} + \text{H}]^+$. Streptomycin: 2 $[\text{M} + \text{H}]^{2+}$, 3 $[\text{M} + \text{H}_2\text{O} + \text{H}]^{2+}$, 4 $[\text{M} + \text{H}]^+$, 5 $[\text{M} + \text{H}_2\text{O} + \text{H}]^+$, 6 $[\text{M} + \text{MeOH} + \text{H}]^+$. Background, b. Spray solvent: 1:1 (v/v) methanol/water with 0.1% (v/v) formic acid. Spray voltage: 5 kV

distinguished at sample loadings of 100 ng (Figure 7B) and 50 ng (Figure 7C): $[\text{M} + \text{H}]^{2+}$ (peak 2), $[\text{M} + \text{H}_2\text{O} + \text{H}]^{2+}$ (peak 3), $[\text{M} + \text{H}]^+$ (peak 4), $[\text{M} + \text{H}_2\text{O} + \text{H}]^+$ (peak 5), and $[\text{M} + \text{MeOH} + \text{H}]^+$ (peak 6). The paper-spray mass spectrum (Figure 7A) shows many background ions that arise from the untreated filter paper, unlike what is observed in polymer-spray (Figures 7B and 7C), but it is possible to remove many of these interferences by a cleaning procedure.²⁸

The $[\text{M} + \text{H}]^+$ (m/z 582.27) was chosen as the precursor ion for tandem mass spectrometry of streptomycin. The fragmentation of the precursor ion to form stable product ions at different collision energies is shown in Figure 8. Ring cleavage occurs at the glycosidic oxygen. The observed product ions are similar to what has been reported using liquid chromatography/tandem mass spectrometry.²⁷

3.4 | Quantitative analysis of streptomycin by OSX polymer-spray

Quantitative analysis of streptomycin in urine by polymer-spray MS was achieved using the peak intensity ratio of the $[\text{M} + \text{H}]^+$ (m/z 582.27) to that of the corresponding $[\text{M} + \text{H}]^+$ peak (m/z 584.29) of dihydrostreptomycin added as an internal standard. The amount of streptomycin in urine was quantified by preparing solutions spiked with streptomycin at different concentrations. The calibration curve was linear over a range of 20 $\mu\text{g}/\text{mL}$ to 500 $\mu\text{g}/\text{mL}$ of streptomycin in urine. The regression correlation coefficient (R^2) was 0.998. The analysis of streptomycin in three urine samples is shown in Table 4. Mean recovery of streptomycin from samples 1 ($n = 3$) and 2 ($n = 3$) were 97.5%




FIGURE 8 Tandem mass spectra of streptomycin (50 $\mu\text{g/mL}$) in water at normalized collision energies: A, 15%; B, 17%; C, 25%; and D, 30%. Spray solvent: 1:1 (v/v) methanol/water with 0.1% (v/v) formic acid. Spray voltage: 5 kV

TABLE 4 Quantitation of streptomycin (m/z 582) in urine, using OSX polymer. Spray solvent: 1:1 (v/v) methanol/water with 0.1% (v/v) formic acid

Substrate	Sample 1 (n = 3) $\mu\text{g/mL}$	Sample 2 (n = 3) $\mu\text{g/mL}$	Sample 3 (n = 3) $\mu\text{g/mL}$
OSX	39.0 \pm 4.9	28.5 \pm 1.5	31.3 \pm 1.6
Actual concentration	40	30	25

Spray voltage: 5 kV. Sample loading: 1 μL .

Dihydrostreptomycin as internal standard (m/z 584).

(RSD 12.6%) and 95% (RSD 5.3%), respectively. The mean recovery from sample 3 (n = 3, RSD 5.1%) was determined to be 20% higher than the actual concentration, which is the lowest concentration of the three samples. Of the administered dose of streptomycin for treatment of urinary tract infections, 60–80% of the dose is excreted in the urine within the first 4 h so, for a dose of 2–4 g/mL, 25–5000 $\mu\text{g/mL}$ in urine can be detected.²⁹ A concentration of 50 $\mu\text{g/mL}$ is considered to be toxic.¹

4 | CONCLUSIONS

We have shown that an organosiloxane polymer is effective as a sample substrate and ESI emitter in the analysis of drug samples without pretreatment by ESI-MS. When compared to paper, organosiloxane polymer substrates allow for mass spectra of a hydrophilic drug at low concentrations not detectable with paper-spray MS. Organosiloxane

polymers achieve a more sensitive analysis, allowing for more accurate quantitation of both hydrophobic and hydrophilic drug compounds. The ability to tailor the polarity and porosity of the polymer allows for the synthesis of a wide range of polymers, and thus opens many possibilities for further development and applications.

ACKNOWLEDGEMENTS

We thank A. C. Mody for help with the MS measurements. This project was supported by the National Institutes of Health STTR Award No. 1R41 GM113337-01 and the National Institutes of Health Award No. R01 MH112188.

ORCID

Maria T. Dulay <http://orcid.org/0000-0002-4678-3480>

Richard N. Zare <http://orcid.org/0000-0001-5266-4253>

REFERENCES

- Mitic SS, Miletic GZ, Kostic DA, Rasic ID. A spectrophotometric study of streptomycin effect on the clinical urea determination. *Chin J Chem.* 2011;29:135-142.
- Streptomycin. Drugs.com website. Available: <https://www.drugs.com/pro/streptomycin.html> (accessed June 28, 2017).
- Brown K, Zaytsoff SJM, Uwiera RRE, Inglis GD. Antimicrobial growth promoters modulate host response in mice with a defined intestinal microbiota. *Sci Rep.* 2016;6:1-13.
- Buszewski B, Noga S. Hydrophilic interaction liquid chromatography (HILIC) - A powerful separation technique. *Anal Bioanal Chem.* 2012;402(1):231-247.
- Notari S, Mancone C, Sergi M, et al. Determination of antituberculosis drug concentration in human plasma by MALDI-TOF/TOF. *Life.* 2010;62(5):387-393.
- So PK, Hu B, Yao ZP. Electrospray ionization on solid substrates. *Mass Spectrom (Tokyo).* 2014;3(S0028):1-8.
- Wang H, Liu J, Cooks RG, Ouyang Z. Paper spray for direct analysis of complex mixtures using mass spectrometry. *Angew Chem Int Ed.* 2010;49:877-880.
- Liu J, Wang H, Manicke NE, Lin JM, Cooks RG, Ouyang Z. Development, characterization, and application of paper spray ionization. *Anal Chem.* 2010;82(6):2463-2471.
- Shen L, Zhang J, Yang Q, Manick NE, Ouyang Z. High throughput paper spray mass spectrometry analysis. *Clin Chem Acta.* 2013;420:28-33.
- Manicke NE, Bills BJ, Zhang C. Analysis of biofluids by paper spray MS: Advances and challenges. *Bioanalysis.* 2016;8(6):589-606.
- Chipuk JE, Gelb MH, Brodbelt JS. Surface-enhanced transmission mode desorption electrospray ionization: Increasing the specificity of ambient ionization mass spectrometric analysis. *Anal Chem.* 2010;81(1):16-18.
- Pasilis SP, Kertesz V, Van Berkel GJ. Surface scanning analysis of planar arrays of analytes with desorption electrospray ionization-mass spectrometry. *Anal Chem.* 2007;79(15):5956-5962.
- Maher S, Jjunju FPM, Damon DE, et al. Direct analysis and quantification of metaldehyde in water using reactive paper spray mass spectrometry. *Sci Rep.* 2016;6(35643):1-10.
- Wang Q, Zheng Y, Zhang Z, Han X, Wang T, Zhang Z. A silica coated paper substrate: Development and its application in paper spray mass spectrometry for rapid analysis of pesticides in milk. *Analyst.* 2015;140:8048-8056.
- Wang X, Zheng Y, Wang T, Xiong X, Fang X, Zhang Z. Metal-organic framework coated paper substrates for paper spray mass spectrometry. *Anal Methods.* 2016;8:8004-8014.
- Ren Y, Wang H, Liu J, Zhang Z, McLuckey MN, Ouyang Z. Analysis of biological samples using paper spray mass spectrometry: An investigation of impacts by the substrates, solvents and elution methods. *Chromatographia.* 2013;76(19-20):1339-1346.
- Damon DE, Maher YS, Yin M, et al. 2D wax-printed paper substrates with extended solvent supply capabilities allow enhanced ion signal in paper spray ionization. *Analyst.* 2016;14:3866-3873.
- Han F, Yang Y, Ouyang J, Na N. Direct analysis of in-gel proteins by carbon nanotubes - Modified paper spray ambient mass spectrometry. *Analyst.* 2015;140:710-715.
- Narayanan R, Sarkar D, Cooks RG, Pradeep T. Molecular ionization from carbon nanotube paper. *Angew Chem Int Ed.* 2014;53:5936-5940.
- Zheng Y, Zhang X, Yang H, et al. Facile preparation of paper substrates coated with different materials and their applications in paper spray mass spectrometry. *Anal Methods.* 2015;7:5381-5386.
- Zhang Z, Xu W, Manicke NE, Cooks RG, Ouyang Z. Silica coated paper substrate for paper-spray analysis of therapeutic drugs in dried blood spots. *Anal Chem.* 2012;17:931-938.
- Damon DE, Davis KM, Moreira CR, Capone P, Cruttenden R, Badu-Tawiah AK. Direct biofluid analysis using hydrophobic paper spray mass spectrometry. *Anal Chem.* 2016;88(3):1878-1884.
- Danks AE, Hall SR, Schnepf Z. The evolution of 'sol-gel' chemistry as a technique for materials synthesis. *Mater Horiz.* 2016;3:91-112.
- Drugbank version 5.0 website. Available: <https://www.drugbank.ca> (accessed June 28, 2017).
- Hayase G, Kanamore K, Fukuchi M, Kaji H, Nakanishi K. Facile synthesis of marshmallow-like macroporous gels usable under harsh conditions for the separation of oil and water. *Angew Chem Int Ed.* 2013;52:1986-1989.
- Jiang J, Zhang H, Li M, et al. Droplet spray ionization from a glass microscope slide: Real-time monitoring of ethylene polymerization. *Anal Chem.* 2015;87:8057-8062.
- Alechaga E, Moyano E, Galceran MT. Simultaneous analysis of kasugamycin and streptomycin in vegetables by liquid chromatography-tandem mass spectrometry. *Anal Methods.* 2015;7:3600-3607.
- Su Y, Wang H, Liu J, Wei P, Cooks RG, Zheng O. Quantitative paper spray mass spectrometry analysis of drugs of abuse. *Analyst.* 2013;138(16):4443-4447.
- Hewitt WL. Treatment of urinary tract infections with streptomycin. *Am J Med.* 1947;2(5):474-484.

SUPPORTING INFORMATION

Additional Supporting Information may be found online in the supporting information tab for this article.

How to cite this article: Dulay MT, Zare RN. Polymer-spray mass spectrometric detection and quantitation of hydrophilic compounds and some narcotics. *Rapid Commun Mass Spectrom.* 2017;31:1651-1658. <https://doi.org/10.1002/rcm.7952>